



CosmoLock® Pedicle Screw System MIS

The CosmoLock® Pedicle Screw System is designed to enhance the surgeon's minimally invasive experience while treating complex pathologies. The system features ergonomic, sleek and intuitive instrumentation focused on reducing common intraoperative challenges encountered by physicians.

CosmoLock® is a registered trademark of the USPTO.



FEATURES

Extremely Low Profile
Breakaway Tower
Tabs Enabling Smaller
Incisions



Bullethead
Rods for Easy
Navigation
Through Soft
Tissue

Dual Lead Threads
With Tapered Tip for
Faster Insertion

4.75, 5.5, 6.5, 7.5, &
8.5mm Self-Tapping
Screws



- Curved Rods Available 35mm to 130mm
- Straight Rods Available 130mm to 250mm in 20mm Increments



Low Profile
Ergonomic Rod
Inserters - Fixed
or Articulating



Innovative Locking
Cap Starter with
Free Spinning Ring
Prevents Head Splay
and Cross-Threading

15mm
Additional
Reduction

Up to 70°
Angulation